

How to Run the model

To run and model and re-produce the best results thins steps need to perform.

Step

1. Run Notebook “[MISA Project PreProcesing Step\(1\) Registration.ipynb](#)” to perform the registration of the Volumes to MNI template.
2. Run Notebook “[MISA Project PreProcesing Step\(2\) Normalization.ipynb](#)” to Perform Preprocessing (Pre-processing pipeline-2 mentioned in report) and to create the excel files that containing the path of the training , validation and testing data. Network Read the data from excel files that have the path of the data.
3. **Folder “Model” Contain the pretrained model and weights.**
4. To run the code please the command “**python train.py --config config_spm_tissue.json**”
5. In the file “**config_spm_tissue.json**” to maintain and configure model

model_path": put ur model weights path (spm_tissue folder)

6. To prepare the Testing Data and After segmentation to bring it back to the original spacing use this Notebook “[PreparingTestingData.ipynb](#)”
7. To run the testing ““**python deploy.py --config config_spm_tissue.json**””
8. Finally to Compute the Dice and Box plot Run the “**Evaluation_MISA_Project.ipynb**”

Resources:

All codes are available in Github: <https://github.com/fitushar/Brain-Tissue-Segmentation-Using-Deep-Learning-Pipeline-NeuroNet>

Best Trained Model Weights Drive Link:

<https://drive.google.com/file/d/1v8gbgLNk5ekgZ5OUgRPYgmN5TKFiuIEU/view?usp=sharing>